AngioVac Circuit

Directions For Use	3
Instrucciones de uso	7
Mode d'emploi	11
Gebrauchsanweisung	15
Istruzioni per l'uso	19
Gebruiksaanwijzing	23
Brugsanvisning	27

16600503-01

D

2023-05

AngioVac Circuit

ENGLISH

Rx ONLY

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

WARNING

Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found, call your sales representative. Inspect prior to use to verify that no damage has occurred during shipping.

For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.

DEVICE DESCRIPTION

The *AngioVac** Circuit is designed to be used with the AngioVac Cannula. It is not designed for long term use (greater than six hours). The circuit includes tubing, connectors, and clamps. The AngioVac Circuit is provided in multiple packaging configurations:

- Circuit with Angio-Sac waste bag
- Circuit with Angio-Sac waste bag and bubble trap

PRIME VOLUME: The AngioVac Circuit has a priming volume of approximately 610 mL

INDICATIONS FOR USE

This circuit is indicated for use in procedures requiring extracorporeal circulatory support for periods of up to six hours.

INTENDED USE

The AngioVac Circuit is intended to be used with the AngioVac Cannula and commonly available extracorporeal equipment (e.g., centrifugal pump, bubble trap (filter) and reinfusion cannula).

CONTRAINDICATIONS

Refer to the AngioVac Cannula Directions for Use (DFU) for procedure-specific contraindications.

WARNINGS

- Inspect product prior to, during, and after use, to ensure that no damage to the product has occurred. Failure to inspect could result in injury to patient or user.
- Carefully examine the circuit for leaks before and during use. Leakage may result in loss of sterility, blood loss or air embolism. If leakage is observed replace the circuit, the leaking component or tighten the leaking connection.
- Instructions for use and manuals for all related extracorporeal circulatory devices should be read prior to use, including but not limited to the AngioVac Cannula, centrifugal pump head, control console and bubble trap.
- As with all medical devices, this device is to be used by or under the direction of trained physicians only.

- User assumes responsibility for any variations in the extracorporeal circulation/procedures that could compromise the intended use of this device.

PRECAUTIONS

- The AngioVac Circuit can only be used in conjunction with the AngioVac Cannula.
- Initiate adequate systemic anticoagulation therapy prior to utilization of this product for the conduct of extracorporeal circulation. A strict anticoagulation protocol should be followed and anticoagulation should be carefully monitored during the procedure.
- Check all pumps and components to ensure that they are set up for proper flow direction prior to initiating extracorporeal circulation.
- Do not allow organic solvents, anesthetic agents or alcohol or alcohol-based fluids to come in contact with the circuit components as they may compromise structural integrity.
- All connections should be finger-tight only. Over-tightening connections may result in cracks of the components and leaks in the circuit.
- To prevent potential for air embolism, always ensure adequate priming and complete de-airing of the circuit prior to use.
- When infusing solutions from bags, remove all air from the bag during setup to prevent air from entering the circuit.
- Carefully monitor circuit for both inflow and outflow obstruction/ occlusion during use.
- During the performance of extracorporeal circulation, do not exceed appropriate pressure ratings of the circuit (520mmHg/10psi).

ADVERSE EVENTS

Possible complications include those normally associated with anticoagulation and extracorporeal circulation.

This product, as with all extracorporeal circulatory systems, has possible side effects, which include but are not limited to infections, blood loss/blood trauma, thrombus formation and embolic events. These may occur if the instructions for use are not followed.

Refer to AngioVac Cannula DFU for additional listing of potential adverse events.

HOW SUPPLIED

Contents supplied STERILE using a Ethylene Oxide (EO) process. Store in a cool, dry, place. Do not use if package is opened or damaged. Do not use if labeling is incomplete or illegible.

OPERATIONAL INSTRUCTIONS

1. The patient should be prepared and draped in the usual and typical sterile manner for percutaneous/surgical vascular procedures.
2. Sterile vascular cannulation should be performed utilizing typical percutaneous or open surgical techniques. Carefully inspect the packaging for any damage or an opened seal.
3. Open barrier and using sterile technique deliver contents to sterile back table.
4. Visually inspect the circuit for any obvious defects in materials or assembly and ensure all protective caps are in place.
5. Gently connect “Y” circuit hub to female connector until it clicks.
6. Remove reinfusion cannula connection and spare parts pouch and place on table.
7. Fully insert tuohy adaptor into red tuohy.
8. Tighten red tuohy until finger tight.
9. Tighten tuohy adaptor until finger tight.
10. Bring contents to patient drape area.
11. Remove tubing from top tray and discard top tray.
12. Remove tubing from bottom tray. Ensure all components have been removed and discard bottom tray.
13. Transfer loose ends of circuit tubing and waste line “T” connection to perfusionist.
14. Remove blue vented cap from circuit tubing and connect to inlet barb connection of bubble trap.

15. Remove white vented cap from waste line "T" connection and connect to outlet barb of bubble trap.
16. Remove remaining white vented cap from waste line "T" connection and connect to inlet barb of the pump head ensuring outlet of pump head is aligned with inlet of bubble trap.
17. Remove red vented cap from circuit and connect to the outlet barb of the pump head.
18. Close waste line tubing clamp.
19. Close the priming line tubing clamp and spike priming fluid bag.
20. Place clamp on blue marker at bubble trap inlet.
21. To fully prime the circuit perform the following:
 - a) Place clamp on blue circuit tubing marker.
 - b) Open priming line tubing clamp to initiate priming.
 - c) Apply pressure to priming fluid bag.
 - d) Slowly open stopcock on bubble trap.
 - e) Continue to apply pressure until fluid has filled circuit tubing, pump head and bubble trap.
 - f) Open waste line clamp to prime waste line and close waste line clamp.
 - g) Close priming line tubing clamp.
 - h) Close stopcock on bubble trap.
 - i) Release clamp off blue circuit tubing marker.
22. Connect pump head to pump per manufacturer's instructions.
23. Start centrifugal pump.
24. Visually inspect entire circuit for leaks.
25. Ensure all air migrates to bubble trap.
26. Shut off pump and open stopcock on bubble trap to remove any remaining air and close stopcock.
27. If desired, connect Angio-Sac waste bag to waste line assembly by removing the waste line end cap and connecting the male luer on the waste bag to the female luer of the waste line until finger tight.
28. Remove red vented cap from the short reinfusion cannula connector and connect to barb of previously inserted reinfusion cannula.
29. Using perfusionist clamps, clamp circuit tubing on both sides of hub assembly on the red and blue markers.
30. Disconnect "Y" circuit hub from reinfusion end of circuit.
31. Connect reinfusion end of circuit (red marker) to reinfusion cannula connector ensuring no air is introduced into the circuit.
32. Remove clamp from reinfusion side of circuit (red marker), operate priming fluid clamp to flush reinfusion cannula with quantity sufficient to clear blood from reinfusion cannula.

WARRANTY

Angiodynamics warrants that reasonable care has been used in the design and manufacture of this instrument. **This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose.** Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond Angiodynamics control directly affect the instrument and the results obtained from its use. Angiodynamics's obligation under this warranty is limited to the repair or replacement of this instrument and Angiodynamics shall not be liable for any incidental or consequential loss, damage or expense directly or indirectly arising from the use of this instrument. Angiodynamics neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. **Angiodynamics assumes no liability with respect to instruments reused, reprocessed, resterilized, modified or altered in any way and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.**

* Angiodynamics, the Angiodynamics logo and AngioVac are trademarks and/or registered trademarks of Angiodynamics, Inc., an affiliate or a subsidiary. All other trademarks are property of their respective owners.

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© 2023 AngioDynamics, Inc. or its affiliates. All rights reserved.
AngioDynamics, Inc.