



Attention - Before using this device, Consult Instructions for Use

## PHOENIX GUIDEWIRE

### INSTRUCTIONS FOR USE



#### DEVICE NAMES

Phoenix Guidewire, Floppy Tip, Light Support

**CAUTION: FEDERAL (US) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON ORDER OF A PHYSICIAN.**

**CAREFULLY READ ALL INSTRUCTIONS PRIOR TO USE, NOTING ALL WARNINGS AND PRECAUTIONS. FAILURE TO DO SO MAY RESULT IN COMPLICATIONS.**

#### INDICATIONS

The Phoenix Guidewire family is designed for percutaneous entry into the peripheral vascular system as an intravascular device by gently negotiating the vascular system while maintaining enough strength and rigidity to enable a catheter to be directed over the guidewire with ease.

#### CONTRAINDICATIONS

This device is not intended for use in the coronary or neurovascular arteries.

#### DESCRIPTION

The Phoenix family is an 0.014 inch diameter X 300 cm length guidewires made of a Silicone coated Nitinol core wire that is tapered at the distal tip where a coil is secured to the distal end. The distal coil is 5.5 cm and consist of Palladium. The tip and wire support of the guidewire are found on product labels and are summarized in Table 1.

**Table 1: Phoenix Guidewire Models and Configurations.**

Catalog No.	Tip and Wire Support Configurations	Wire Tip	Core Material	Diameter	Length
PG14300LF	Soft Tip /Light Support	5.5 cm	Nitinol	0.36mm (0.014")	300 cm

#### WARNINGS

- The Phoenix Guidewire is supplied sterile and intended for one-time use only. Do NOT resterilize and/or reuse. Resterilization or reuse may potentially compromise device performance and safety and may increase the risk of infection.
- Use the Phoenix Guidewire products prior to the "Use By" ("Expires") date specified on each individual package label.
- When the Phoenix Guidewire is exposed to the vascular system, it should not be advanced or retracted except under direct fluoroscopic observation. If resistance is met during advancement or retraction, determine the cause of the resistance before continuing.
- Do not withdraw guidewires through a metal trocar or metal needle.

- Do not advance or withdraw guidewires against resistance until the cause of the resistance has been determined. Excessive force against resistance may result in damage to guidewire and catheter or vessel perforation.
- Volcano does not recommend a particular technique for the use of this device. The steps contained in the preceding directions discuss the Seldinger Technique of percutaneous entry are for information purposes only. Each physician should evaluate their appropriateness according to individual patient condition and his or her medical training and experience.
- Phoenix guidewires incorporate various performance characteristics from tip shape to body stiffness. Extreme care must be taken in choosing the correct guidewire or to minimize the risk of vessel perforation or vascular damage.
- Guidewires by nature of their construction will collect blood and other foreign matter in the coiled tip. No type of cleaning will completely remove this material. Therefore, they are intended for single patient or one time use only. DO NOT RESTERILIZE.
- Do not use, or attempt to correct, a Phoenix Guidewire if it is bent or kinked or has any evidence of damage as this may result in breakage or compromised performance

#### CAUTIONS

- Prior to use, the Phoenix Guidewire should be inspected damage. If damaged, DO NOT USE.
- For use only by physicians trained in percutaneous peripheral interventional procedures.
- Carefully read all instructions prior to use. Observe all warnings and cautions. Failure to do so may result in complications

#### ADVERSE EFFECTS

Potential adverse effects associated with the use of this device and other interventional catheters include but are not limited to the following (alphabetical order):

- Access site injury, including pain
- Arterial dissection
- Arterial perforation or pseudoaneurysm
- Arterial spasm or abrupt or sub-acute closure
- Death
- Embolism, including thrombus, plaque, air, device, etc.
- Fracture of any component of the device that may or may not lead to serious injury or surgical intervention
- Hematoma requiring surgery
- Myocardial infarction
- Arrhythmia
- Renal Failure
- Myoglobinuria
- Reaction to contrast media, procedure medications or wire materials, including allergic reaction
- Revascularization or vascular injury which may require surgical repair or emergency surgery

The occurrence of these adverse effects may lead to the need for repeat catheterization/angioplasty, emergency bypass surgery, or death.

#### HOW SUPPLIED, STERILIZATION, AND EXPIRATION

The Phoenix Guidewire is packaged and sterilized individually. **NOTE:** Guidewire, length, diameter, tip and wire support configurations are indicated on the product label. The Phoenix Guidewire is sterilized using ethylene oxide gas. Each device is sterile if the sealed pouch is unopened and undamaged. All are intended for single use only and should not be reused or re-sterilized. Use before the "Use by" dated printed on the package labels. Store in a dry, cool place.

#### INSTRUCTIONS FOR PROCEDURE

##### INSPECTION OF PHOENIX GUIDEWIRE COMPONENTS












Prior to use, all equipment to be used for the procedure should be examined carefully for defects. Examine the packaging for cuts, tears, or other breach of the sterile barrier. Examine the Phoenix Guidewire for bends, kinks or other damage. Do not use any defective guidewire



##### DIRECTIONS FOR USE:

1. Inspect guidewire prior to use for tip shape, bends, kinks, coil separation. If damaged, DO NOT USE.
2. Using sterile technique, localize and puncture the vessel with a needle cannula.
3. Remove the needle, leaving the cannula in place.
4. Insert distal end of the guidewire through the cannula and into vessel.
5. Remove the cannula, leaving the guidewire within the lumen of the vessel.
6. Pass the catheter over the guidewire within the lumen of the vessel under fluoroscope guidance to the desired position.
7. With the guidewire in place, follow the instructions for use provided by the manufacture of the introducer, catheter or other devices, as well as the medical facility's standard procedure for these types of devices.

##### DEVICE DISPOSAL:

The Phoenix Guidewire is disposable per hospital biohazard procedures.

 Use Before Date	 Nonpyrogenic
 Store in a dry, dark, cool place	 Content: One (1)
 Do not use open or damaged packages	 Do not Re-Sterilize
 Catalogue number	 Batch Code
 Sterilized using Ethylene Oxide	 Single Use Only
 Prescription Only	

 Manufactured for:	Volcano AtheroMed, Inc. 1455 Adams Dr. #1120 Menlo Park CA 94025 USA (800) 228-4728 (916) 638 8112 fax	 Authorized Representative	Volcano Europe BVBA/SPRL Excelsiorlaan 41 B-1930 Zaventem, Belgium +32.2.679.1076 +32.2.679.1079 fax
--	---	---	---

©2016 Volcano AtheroMed, Inc. All rights reserved.

ATHEROMED, PHOENIX, and the ATHEROMED and PHOENIX logos are trademarks or registered trademarks of Volcano AtheroMed, Inc. in the U.S. and other countries.

The PHOENIX ATHERECTOMY System and PHOENIX GUIDEWIRE (including components and/or methods thereof) may be protected under one or more pending United States Patent Applications, as well as corresponding international patent applications.

#### LIMITED WARRANTY AND DISCLAIMER

Volcano AtheroMed, Inc. warrants that reasonable care has been used in the design and manufacture of this medical device.

Handling, storage and preparation of this medical device as well as other factors relating to the patient, diagnosis, treatment, surgical procedures and other matters beyond the control of Volcano AtheroMed, Inc. directly affect this medical device and the results obtained from its use. Further, no representation or warranty is made that a Volcano AtheroMed, Inc. product will not fail. Volcano AtheroMed, Inc. disclaims responsibility for any medical complications – including death – directly or indirectly resulting from the use of this product. Except as expressly provided by this limited warranty, VOLCANO ATHEROMED, INC. IS NOT RESPONSIBLE FOR ANY DIRECT, INCIDENTAL OR CONSEQUENTIAL DAMAGES BASED ON ANY DEFECT, FAILURE OR MALFUCTION OF ITS PRODUCTS WHETHER THE CLAIM IS BASED ON WARRANTY, CONTRACT, TORT OR OTHERWISE.

Volcano AtheroMed, Inc. assume no liability with respect to medical devices reused, reprocessed or resterilized and makes no warrants – expressed or implied – including but not limited to merchantability of fitness for a particular purpose with respect to such medical device.

THIS WARRANTY IS IN LIEU OF AND EXCLUDES ALL OTHER WARRANTIES NOT EXPRESSLY SET FORTH HEREIN, AND VOLCANO ATHEROMED, INC. MAKES NO WARRANTY – EXPRESS OR IMPLIED – INCLUDING BUT NOT LIMITED TO ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE.